

arkray

FOR QUICK & EASY HEMOGLOBIN TESTING

PocketChem™ Hemo G

Handheld Point-of-Care Hemoglobin Meter



POC Reports App



Scan QR Code to Download
(Android)

www.arkray.co.jp/english

www.arkray.co.in

● SMALL

- Small sample Size - 7 μ l
- No large storage space required
- Most of the area covered by LCD for better view

● SIMPLE

- Easy 3 Step Process
- Results in 5 seconds
- Whole blood / Venous sample
- Special Micro-pipette to collect exactly 7 μ l sample
- LED illumination to mark blood application area

● SPECIFIC

- Designed for Hb & Hct measurement only
- Accuracy calibrated to cell counter

● MICRO-PIPETTE

- Designed to collect 7 μ l blood only
- Disposable Micro-pipette tips
- Each bottle of strips provided with equal no. of Micro-pipette tips



EASY 3 STEP PROCESS



1

Insert strip to turn on the meter



2

Collect 7 μ l sample using special micro-pipette



3


Apply sample to the application area on the strip



Read results in 5 seconds



For more information contact:



Customer Support: +91 77770 53311
(Mon - Fri 9.30am to 6.30pm)
Email: customersupport@arkray.co.in

ARKRAY Healthcare Pvt. Ltd. 701-702, Opulence, 6th Road, TPS III, Santacruz (E), Mumbai - 400055, Maharashtra, India.