

FOR TRUSTWORTHY BP MONITORING AT HOME

TRUSTCHECKTM BPM^{2.0}

DIGITAL UPPER ARM BLOOD PRESSURE MONITOR

5 YEARS^{*}
EXTENDED
WARRANTY



- WHO INDICATOR
- UNIVERSAL CUFF
- QUICK INFLATION
- 120 TESTS MEMORY FOR 2 USERS
- IRREGULAR HEARTBEAT DETECTION
- FACILITY FOR MICRO USB POWER ADAPTER

M.R.P: Rs.2045 incl. of all taxes

TRUSTCHECKTM BPM^{2.0}

DIGITAL UPPER ARM BLOOD PRESSURE MONITOR

MODEL NO.: B56



WHO INDICATOR

Makes it easy to understand results



UNIVERSAL CUFF

Fits multiple arm sizes



QUICK INFLATION

For faster testing time



120 TEST MEMORY FOR 3 USERS

to store readings of every user separately



FACILITY OF USB POWER ADAPTER

To use the meter without batteries



TALKING METER

To make testing process easier



IRREGULAR HEARTBEAT DETECTION

Can help in detection of pulse irregularities at an early stage.

TECHNICAL SPECIFICATION

Product Type	Upper arm type-B56
Measurement Range	DIA: 40-130mmHg; SYS: 60-230mmHg Pulse 40-199beats/min
Accuracy	Blood pressure: ± 3 mmHg Pulse: $\pm 5\%$ of the reading
Memory	2*120 groups memory (Double users)
Average Function	Last 3 groups average measuring value
Material	ABS+LCD display
Cuff Circumference	22-40cm
Power Source	Internal-DC 6V (4*AAA) / External-DC 5V 1A Temperature: 5°C to 40°C
Operating Conditions	Humidity: 15% to 93%RH Pressure altitude: 70Kpa~106Kpa Temperature: -25°C to 70°C
Storage And Shipping Conditions	Humidity: $\leq 93\%$ RH Pressure altitude: 70Kpa~106Kpa

UNIVERSAL CUFF

TRUSTCHECK BPM 2.0
FITS 20 X 40 cm arm size



OTHER BRANDS
FIT 20 X 32 cm arm size

For more information contact:

 +91 777705331
customersupport@arkray.co.in